

# PerFect® TCS II

The state-of-the-art electrosurgical system for precise and efficient tissue contouring


Solid state circuitry eliminates warm-up time

Waterproof footswitch for hands-free electrode activation

Single control for easy mode selection and intensity adjustment

Swivel cord connector improves handling

Dispersive electrode completes energy flow path



8 pcs. Autoclavable electrode sheaths

## Ordering Information

PerFect® TCS II	230V
Cat. No.	S8230CE
	S8230UK


Included Components	Qty.
PerFect TCS II Unit 230V	1
Cord Assembly S6008CE	1
Electrode Set S7010A (S6011A, S6012A, S6015A)	3
Dispersive Electrode S213	1
Footswitch assembly S5006	1
Handpiece Holder S7001	1
Introductory DVD	1
Step-by-step Owner's Guide	1

Cat. No.	Electrode Sheaths	Qty.
S6010A	Complete Set S6011A thru S6018A	8
S7010A	Included Set S6011A, S6012A, S6015A	3
S6011A	Coag Ball	2
S6012A	45° Straight Knife	2
S6013A	45° Tapered Knife	2
S6014A	Straight Round Loop	2
S6015A	45° Long Loop	2
S6016A	Straight Long Loop	2
S6017A	45° Diamond Loop	2
S6018A	90° AP 1½	2

## Technical Data

Measure: 108 × 159 × 222 mm  
Weight: 2.4 kg

© 2013 Coltène/Whaledent AG – www.coltene.com



**Coltène/Whaledent AG**  
Feldwiesenstr. 20  
9450 Altstätten/Switzerland  
Tel. + 41 71 757 53 00  
Fax + 41 71 757 53 01  
info.ch@coltene.com

**Coltène/Whaledent Inc.**  
235 Ascot Parkway  
Cuyahoga Falls,  
Ohio 44223/USA  
Tel. + 1 330 916 8800  
Fax + 1 330 916 7077  
info.us@coltene.com

**Coltène/Whaledent Ltd.**  
The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF / UK  
Tel +44 1444 235486  
Fax +44 1444 870640  
info.uk@coltene.com

 **COLTENE®**